

Experience the Depth
Created by Two Frequencies
The Intense Harmony of
6.78 MHz and 2 MHz

XERF combines these two frequencies to deliver thermal energy that penetrates deeply, transmitting heat to more skin layers without the need for invasive procedures.



FIT MY SKIN, FIT MY MOMENT

Customized Treatment Just for You with Precise Depth and Level

Experience XERF's fine-tuned 3 levels of Depth and 10 levels of Intensity. Each customized treatment delivers heat at varying depths, benefiting your skin in a way tailored to its characteristics.

Accurate Impedance Feedback

Real-time impedance feedback analyzes your skin's condition before the procedure, enabling a treatment tailored to fit your skin.

How long does the treatment take?

XERF can be performed immediately without any anesthesia. Depending on the number of shots, it takes 15 to 30 minutes.

Can I resume my daily activities right after the treatment?

XERF is a zero-downtime treatment that delivers energy deep into the skin. You can generally resume your daily activities immediately after the treatment without any recovery period.

Will it hurt?

XERF's ICD Cooling effectively reduces pain by lowering the epidermal surface temperature. The Wave Fit Pulse integrated into XERF also ensures a stable treatment.

Is it safe?

XERF is an approved medical device by the Korean MFDS (Ministry of Food and Drug Safety), which confirms its safety.

Which XERF EFFECTOR and how many treatments will I need?

The appropriate number of shots and treatment sessions vary depending on your individual skin condition and treatment area. Please consult a nearby clinic for your personalized XERF treatment plan.

XERF

Meet the New
Monopolar RF : XERF

SHALLOW
MIDDLE
DEEP



DEPTH

CYNO|SURE
LUTRONIC

Advertisement Review Number: Universal Electrosurgical Device Association
-2024-33-051 (Validity Period: 2027-08-27) This product is a medical device.
Please read and follow the precautions and instructions for use carefully.
Copyright©Lutronic Corporation. All rights reserved.

XERF ambassador
Shinhye Park

FIT THE

WAVE FIT

FIT WITH TRUST

Safety through
Technology

Constantly Evolving
Wave Fit Pulse


XERF's Wave Fit Pulse adapts the RF pulse in alignment with the amount of energy delivered, providing stable treatments with the perfect balance of energy and cooling effect.

PULSE


Diverse Line-up

We offer a variety of EFFECTOR sizes to fit a broad range of treatment areas. Choose the EFFECTOR that suits your needs from the 4 options available.

EFFECTOR 60



EFFECTOR 40



EFFECTOR 10



EFFECTOR 05



EFFECTOR Surface Temperature Monitoring System

The surface temperature of the treatment area is monitored in real-time to prevent overheating. If the surface temperature exceeds 43°C, the RF output is automatically shut off.

Stabilizing ICD Cooling

XERF's 3-level ICD Cooling system, tailored to match the treatment intensity, prevents excessive epidermal temperature rise. The perfect temperature balance reduces skin irritation and maximizes comfort.

Fit to Shape
A Presence Felt Anywhere

The XERF EFFECTOR package is inspired by the multi-dimensional inner selves of individuals who always follow their own paths.



Memo >